NAUTIAGENE

Product Information

NautiaZ Plasmid DNA Extraction Mini Kit

(100/300 prep)

Cat. No. : NPZ-S100/NPZ-S300 Sample : 1-4 ml of bacterial cultures Yield : Up to 50 µg

NPZ-S100	NautiaZ Plasmid DNA Extraction Mini Kit (100 prep)
NPZ-S300	NautiaZ Plasmid DNA Extraction Mini Kit (300 prep)
NPZ-M020	NautiaZ Plasmid DNA Extraction Midi Kit (20 prep)
NPZ-L010	NautiaZ Plasmid DNA Extraction Maxi Kit (10 prep)
NPZ-L020	NautiaZ Plasmid DNA Extraction Maxi Kit (20 prep)

Contents

	NPZ-S100T	NPZ-S100	NPZ-S300
PS1 Buffer*	1 ml	25 ml	65 ml
PS2 Buffer	1 ml	25 ml	65 ml
PS3 Buffer	1.5 ml	35 ml	95 ml
W1 Buffer	2 ml	45 ml	125 ml
W2 Buffer**	300 µl X2	15 ml	25 ml x2
Elution Buffer	1 ml	10 ml	30 ml
RNase A (50mg/ml)	Added	50 µl	150 ul
PS Column	4 pcs	100 pcs	300 pcs
Collection Tube	4 pcs	100 pcs	300 pcs
User Manual	1	1	1

*Add provided RNase A to PS1 Buffer and store at 4°C. **Add 60 ml / 100 ml ethanol (96-100%) to W2 Buffer prior to the initial use.

Buffer Preparation

• Add ethanol (96-100%) to the Wash Solution prior to first use

	NPZ-S100T	NPZ-S100	NPZ-S300
PS1 Buffer	-	25 ml	60 ml
RNase A		50 μl	150 μl
W2 Buffer	300 μl X2	15 ml	25 ml
ethanol (96 ~ 100%)	1.2 ml X2	60ml	100ml

Important Notes

- 1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
- 2. Check PS2 Buffer before use. If any precipitation formed, warm the buffer in a 37°C waterbath to dissolve.
- 3. To avoid acidification of PS2 Buffer from CO_2 in the air, close the bottle immediately after use.
- 4. All centrifuge steps are done at full speed (14,000 x g) in a microcentrifuge.

Growth of Bacterial Cultures

- 1. Pick a single colony from a freshly streaked selective plate to inoculate 1-4 mL of LB medium supplemented with the appropriate selection antibiotic. Incubate for 16-18 hours at 37°C while shaking at 200-250 rpm. Use a tube or flask with a volume of at least 4 times the culture volume.
- 2. Harvest the bacterial culture by centrifugation at 8000 rpm ($6800 \times g$) in a microcentrifuge for 2 min at room temperature. Decant the supernatant and remove all remaining medium.

Storage

RNase A should be store at 4°C , the other Buffer and Columns store at room temperature.



Transfer the supernatant carefully to PS Column.

STEP PROCEDURE 5 DNA 1 Transfer 1-4 ml of well-grown bacteria culture to a Bindir Harvestmicrocentrifuge tube(not provided). Descend the bacteria by centrifuging at 14,000 x g for 1 ing min and discard the supernatant completely. Add 200 µl PS1 Buffer (RNase A added) to the pellet and 2 Resuspenresuspend the cells completely by pipetting or vortexing. Wash sion Note: Make sure that RNase A has been added into PS1 Buffer whenfirst use. • No cell pellet should be visible after resuspension of the cells. Add 200 µl PS2 Buffer and mix gently by inverting the 3 tube 10 times to lyse the cells and incubate at room Lvsis temperature for 2 mins until the lysate is homologous. Note • Do not vortex, vortex may shear genomic DNA. • Do not proceed this step over 5 min. Add 300 µl PS3 Buffer and invert the tube 10 times 7 immediately but gently. Dry colum Note • Do not vortex, vortex may shear genomic DNA. • Invert immediately after adding PS3 Buffer will avoid asymmetric

Centrifuge at 14,000 x g for 3 mins. During centrifuging, place a PS Column in a Collection Tube.

8 DNA

Elution

g	Place the PS column back into the Collection Tube. Note • Do not transfer any white pellet into the column.
	Add 400 μl W1 Buffer to PS Column. Centrifuge at 14,000 x g for 30 seconds then discard the flow-through. Place the PS column back into the Collection Tube.
	 Add 600 μl W2 Buffer to PS Column. Centrifuge at 14,000 x g for 30 seconds then discard the flow-through. Place the PS column back into the Collection Tube. Note Make sure that ethanol (96-100 %) has been added into Wash uffer when first open.
n	Centrifuge at 14,000 x g for an additional 2 min to dry the column. Note • Important step ! This step will remove the residual liquid completely that will inhibit subsequent enzymatic reaction.
	Place PS Column to a new 1.5 ml microcentrifuge tube (not provided).

Add 50 - 200 µl of Elution Buffer or ddH₂O to the membrane center of PS Column. Stand the column for 2 mins. Note: For effective elution, make sure that the elution solution is dispensed on the membrane center and is absorbed completely. If plasmid DNA is larger than 10 kb, use preheated 70°C Elution Buffer to improve the elution efficiency.

Centrifuge at 14,000 x g for 2 mins to elute plasmid DNA.

9 Pure DNA Store plasmid DNA at 4 °C or -20 °C.

4

Neutralization

precipitation.